



FOX IV ⁸¹⁰ DIODE LASER

A.R.C.
LASER

MADE IN GERMANY ○ ● ● ●

ALL-ROUND DIODE LASER FOR OPHTHALMOLOGY

SMART
MODERN
PORTABLE

GLAUCOMA
SPECIALIST







FOX IV 810

YOUR BEST OPTION FOR



Glaucoma IOP Reduction



Retinal Photocoagulation



Lacrimal Surgery



ENERGY & TIME
DISPLAY MODE




AUDIO IN
5 LANGUAGES




METRONOME
FUNCTION

The ultimate
solution for
CPC & μCPC



The ultimate diode laser for open-angle glaucoma

- **Time or energy counter?**
FOX IV 810 offers both options.
- **Speaking companion**
FOX IV 810 provides a 5 language voice output and counts applied energy or time in steps of 10.
- **Always the right speed!**
FOX IV 810 improves μ CPC use with its metronome function offering an audio signal in 1 s or 1 J intervals.



Next generation portable diode Laser

- Durable metal housing (17.4 x 14.2 x 16.3 cm)
- Large user-friendly colour touch display
- Highly portable, battery powered FOX IV 810 with all accessories in a compact carrying case
- Lightweight - only 1.53 kg

Powerful and reliable

- 810 nm diode laser with power output up to 8 Watts
- Safety-feature: power limitation of 3 W for anterior and posterior treatments
- Color coded application user interface
- Important treatment parameters prominently displayed

Reputation you can trust

- Maintenance-free technology
- Based on smart features of the FOX 810
- More than 15 years of diode laser expertise and thousands of units sold worldwide

COMPLETE RANGE OF ACCESSORIES

Glaucoma

- Glaucoma Probes



Retinal Coagulation

- LIO
- Endo Probes



Lacrimal Surgery

- DCR Probes



www.ARCLASER.COM



CE 0123 A.R.C. Laser ist zertifiziert nach EN ISO 13485:2012 + AG:2012 93/42 EWG

A.R.C.
LASER

MADE IN GERMANY ○ ● ●



Publisher and Copyright

A.R.C. Laser GmbH
Bessemerstr. 14
90411 Nürnberg

Telefon: +49 (0) 911-21779-0
Telefax: +49 (0) 911-21779-99
E-Mail: info@arclaser.de